

Designed to make your life in the laboratory easier!

All NEW
GFD®
Lab Filter
Dryer



The new **GFD**® Bench Top Nutsche
Filter Dryer efficiently performs
laboratory process development



PSL has offices around the globe to support our worldwide customers
+44 (0)151 448 7700 info@powdersystems.com www.powdersystems.com

GFD®

from  psl

Easy laboratory synthesis with direct scale up capability up to production



DATA RECORDING & EXPORT

- HMI Display
- Real Time Values and Graphs
- Export CSV files via USB or Ethernet

PROCESS CONTROL

- HMI/PLC System
- Accurate Process Control
- Sequence Recipe Import
- Full Automation & Repeatability

BENEFITS

- Efficient Filtration & Drying
- Removable Filter Basket
- Easy to Install and Use – Plug and Play
- Enhanced Ergonomics
- Easy Cleaning

OPTIONS

- cGMP Documentation
- Pyrometer and Pressure Relief & Sensor
- Alternative product contact materials
- Range of filtration media and pore sizes
- TCU and/or Vacuum Pump Packages
- Metallic vessel up to 4barg pressure rated
- High Containment solution, GFD FilterBox

Interchangeable sizes and materials of vessel available

	MINI LAB	LAB	MAXI LAB
Vessel Volume	0.3L	2L	10L
Cake Volume	0.1L	0.5L	5L

Borosilicate 3.3 Glass, 316L Stainless Steel or C22 Alloy



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INNOVATION
2016



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INTERNATIONAL TRADE
2016



Serving our customers since 1989